

**New &  
Improving**

*You're  
the Artist*

**SPIDENT**

## Variety of Indications with EsFlow

# EsFlow<sup>®</sup>

## Light Cured Nanohybrid Flowable Composite Resin

### Advantages

- Highly esthetic
- Superior compressive strength
- Low polymerization shrinkage
- Flow on demand handling
- Optimal radiopacity
- Strength equivalent to natural enamel after light-cured
- Available in 2 viscosity : EsFlow<sup>®</sup>, EsFlow<sup>®</sup> LV
- Available shades : EsFlow<sup>®</sup> : A1, A2, A3, A3.5, B1, B2, Blue  
EsFlow<sup>®</sup> LV : A1, A2, A3



SPIDENT CO., LTD.

203 & 312, Korea Industrial Complex, 722, Gojan-Dong, Namdong-Gu, Incheon, Korea 405-821

Tel. 82-32-821-0071 Fax. 82-2-318-0034 E-mail. synergy@spident.co.kr [www.spident.co.kr](http://www.spident.co.kr)

CE  
0120

# EsFlow<sup>®</sup>


## Light Cured Nanohybrid Flowable Composite Resin




### • Indication

- Restoration in Class III~V
- Deep Class I applicable
- Base / Liner under direct restorations
- Undercut blockout
- Usable in Orthodontics
- Repair of composite
- Pit & fissure sealant

### • Physical Properties




### Fluidity




EsFlow (SPIDENT)

EsFlow LV (SPIDENT)



### Slumping



EsFlow (SPIDENT)

Product D

Product E

Product F

### • Delivery Form

Syringe Type

#### EsFlow<sup>®</sup> Syringe Refill (Pack)

Contains : 2g x 2 syringes, 8 Dispensing tips  
Shade : A1, A2, A3, A3.5, B1, B2, Blue



Syringe Type

#### EsFlow<sup>®</sup> LV Syringe Refill (Pack)

Contains : 2g x 2 syringes, 8 Dispensing tips  
Shade : A1, A2, A3

